

# Tip

## Null Modem Cable Configuration Class-PA Server PC to Model 412

PC- SubD 9Pin	Model412 – SubD 9Pin
2	3
3	2
7	8
8	7
5	5

**Guided Wave Incorporated**

3033 Gold Canal Drive  
Rancho Cordova, CA 95670  
USA  
Tel: 916.638.4944; Fax: 916.635.8458  
[gwinfo@guided-wave.com](mailto:gwinfo@guided-wave.com)

[www.guided-wave.com](http://www.guided-wave.com)

**Guided Wave BV**

PO Box 427  
7550 AK Hengelo (o)  
The Netherlands  
Tel: +31.74.2595390; Fax: +31.74.2595752  
[info@guided-wave-europe.com](mailto:info@guided-wave-europe.com)